

# **RX MARINE INTERNATIONAL** Total Solution Total Protection

AN ISO CERTIFIED COMPANY



# **Oxygen Scavenger Plus**

Part/Order no:	Packing
RXSOL-50-5007-20	20 Ltr
RXSOL-50-5007-25	25 Ltr

### **Oxygen Scavenger Plus:**

Oxygen Scavenger Plus is a catalyzed liquid solution of Hydrazine and diethyl hydroxylamine (DEHA). Its volatile properties ensure distribution throughout the boiler and condensate system, and hence protects against oxygen corrosion in all parts of the system. The product also provides the required conditions for the establishment of a passivating layer of magnetite on all inner surfaces. RXSOL-50-5007 is a low toxicity product, suitable for high, medium and low pressure boiler systems.

#### **Product Specifications:**

Density.	1.
Packing,Type.	PLASTIC.
Volume.	25 I.
PH.	10.5.
Color.	LIGHT YELLOW.
Appearance.	LIQUID.

## Directions for use:

This product is formulated for use in conjunction. Oxygen Scavenger Plus can be used in any boiler system, and in conjunction with mechanical desecration systems. Removal of dissolved oxygen is vital for prevention of corrosion and especially pitting corrosion in boilers. When dosed into a boiler system, Oxygen Scavenger Plus will react with dissolved oxygen and form non corrosive compounds. The product is volatile, and left over product from oxygen scavenging in the pre-boiler section will evaporate and assist protection of the steam and condensate system. No solid materials are produced when using Oxygen Scavenger Plus.

# Dosing method:

Rx Marine International recommends continuously into the boiler feed line and Boiler Water Treatment Dosing Unit. Dosing should be controlled so that the DEHAresidual in the condensate is between 0.08-0,30 ppm. The consumption will depend on feed water temperature and the amount of water fed into the boiler. For more information about initial dosage and dose rate please consult your nearest Customer Centre.

Test result condensate PPM DEHA 0.00 - 0.08 0.08 - 0.30 0.30 +Increase dose 25% Maintain dose Decrease dose 25% Sampling and testing.

Arepresentative sample of Condensate should be drawn for analysis daily. The sample should always be taken from the same sampling point, cooled and tested immediately By Rxsol Test Kit. It is important that regular testing is carried out, to ensure the correct level of treatment is maintained

105 - A Wing, BSEL TECH PARK, Plot No. 39/5, Sector-30, Opp. Vashi Railway Station Navi Mumbai 400705 Maharashtra ,India Tel:- (0091) -22- 65113333/5555/ 9999 Email:- mail@rxmarine.com website :- www.rxmarine.com

#### Features:

- Fast acting, liquid oxygen scavenger.
- Safe and easy to use, low toxicity.
- Reduces corrosion of iron and copper, increasing system life and reliability.
- Neutralizes acids occurring in condensate system.
- Volatile product, provides distribution and protection throughout the boiler system, economical in use.
- Suitable for use in conjunction with combined and co-ordinated treatment programmes.
- Organic product, no dissolved solids added.
- Simple test to determine treatment level.

